


standards at a constant atmospheric pressure or negative pressure of 10⁻⁶mBar





Why a vacuum comparator?



After redefinition of the kilogramme unit and based it on a Planck constant, the measurement carried out on the Watt (Kibble) scale was approved by the leading-edge experiment. Due to introducing a new element (measurement on the Watt scale) into the uncertainty budged, the kilogramme is determined with the worse uncertainty than it had place before the redefinition.

Additionally the measurements were to be carried out in the vacuum setting because of the interferometers used. It resulted in the necessity of equipping the national metrology institutes in a new comparator to deal with the national kilogramme standards with the higher resolution of 0,1 ug (in order to decrease the uncertainty of the measurement), equipped also in the vacuum or constant pressure chamber (constant air density).



The extension variants of the constant pressure comparator



AVK-1000.5Y.CP

Automatic Constant Pressure Mass Comparator



AVK-1000.5Y

Automatic Vacuum Mass Comparator equipped with pumps



AVK-1000.5Y.LLS

Automatic Vacuum Mass Comparator equipped with pumps and Load-Lock System



AVK-1000.5Y.CP

Automatic Constant Pressure Mass Comparator

AVK-1000.5Y.CP is a constant-pressure comparator (constant air density) with a maximum load capacity of **1 kg**.

It is supplied with a **6-position mass standard magazine**, a self-centring suspended weighing pan and automatic external adjustment with a (replaceable) standard weight.


The comparator can also be used to compare silica hemispheres with a diameter from 40 to 100 mm, as well as cylindrical standards of ø (22-95)x110.


| Model | Maximum capacity | Readability | Standard repeatability | E0 | E1 | E2 | F1 | F2 |
|--------------------|------------------|-------------|------------------------|--------------|--------------|--------------|--------------|--------------|
| AVK- 1000.5Y.CP | 1002 g | 0,1 μg | 1 µg | 100 g - 1 kg |



Main advantages of the AVK-1000.5Y.CP comparator

6-position mass standard magazine for calibration of cylindrical standards ø (22-95)x110 and silica hemispheres with a diameter from 40 to 100 mm








Main advantages of the AVK-1000.5Y.CP comparator

Self-centring suspended weighing pan

Elimination of a non-centricity error thanks to the use of a special suspended weighing pan.







AVK-1000.5Y

Automatic Vacuum Mass Comparator equipped with pumps

The constant-pressure comparator can be equipped with a first-level vacuum pump and turbomolecular pump as well as a vacuum meter.

The comparator has advantages of the constant-pressure comparator, can generate a vacuum of **10**⁻⁶ **mBar**, and allows performing comparison in the vacuum.

| Model | Maximum capacity | Readability | Standard repeatability | E0 | E1 | E2 | F1 | F2 |
|-----------------|------------------|-------------|------------------------|--------------|--------------|--------------|--------------|--------------|
| AVK- 1000.5Y | 1002 g | 0,1 μg | 0,5 μg | 100 g - 1 kg |



Main advantages of the AVK-1000.5Y comparator

The automatic calibration mechanism with potential replacement of the adjustment weight for a lower-mass weight in order to check sensitivity of the comparator.







Main advantages of the AVK-1000.5Y comparator

It is possible to perform comparison in the atmosphere of inert gases.



Vacuum of up to 10⁻⁷ mBar.

The comparator is equipped with 2 pumps (first-level vacuum pump and turbomolecular pump), as well as a high-class vacuum meter.







AVK-1000.5Y.LLS

Automatic Vacuum Mass Comparator equipped with pumps and Load-Lock System

The vacuum comparator can be additionally equipped with a "LOAD LOCK" system. It is responsible for feeding mass standards and is supplied with a separate pump unit that allows generating a vacuum similar to the one that is present in the main chamber. Thanks to lower volume of the loading chamber, similar vacuum can be generated a way faster (around 4 hours). Afterwards it is possible to load or unload one mass standard from the main chamber.

The **RADWAG's** original "LOAD LOCK" system allows loading the standard from the main chamber into the special vacuum container and transferring the mass standard in the vacuum.

| Model | Maximum capacity | Readability | Standard repeatability | E0 | E1 | E2 | F1 | F2 |
|------------------|------------------|-------------|------------------------|--------------|--------------|--------------|--------------|--------------|
| AVK- 1000.LLS | 1002 g | 0,5 μg | 0,5 μg | 100 g - 1 kg |



Main advantages of the AVK-1000.5Y.LLS comparator

LOAD LOCK system

The comparator is supplied with a special mass standard feeder with a separate pumping unit.

Vacuum containers

The LOAD LOCK system is used to load and unload a mass standard from a dedicated container that allows transferring and storing the mass standard in the vacuum.







Environmental requirements for using the vacuum (constant-pressure) comparator

In view of its enormous reading unit (10 billion units), the vacuum comparator requires specific ambient conditions.

The temperature and vibrations are the main factors influencing the repeatability and final results. Temperature fluctuations: <0,1°C/24h.






Date




Environmental requirements for using the vacuum (constant-pressure) comparator

In view of its enormous reading unit (10 billion units), the vacuum comparator requires specific ambient conditions.

The temperature and vibrations are the main factors influencing the repeatability and final results. Temperature fluctuations: <0,1°C/24h.






Date



AVK 1000.5.CP comparator repeatability results


| Mana | MEASUREMENTS | | | | | | | | |
|--------|--------------|-----------|-----------|-----------|-----------|---------|--|--|--|
| Mass | 1 SD [µg] | 2 SD [µg] | 3 SD [µg] | 4 SD [µg] | 5 SD [µg] | SD [µg] | | | |
| 1000 g | 0,86 | 0,52 | 0,43 | 0,92 | 0,74 | 0,694 | | | |
| 500 g | 0,34 | 0,77 | 0,6 | 0,83 | 0,54 | 0,616 | | | |
| 200 g | 0,63 | 0,32 | 0,29 | 0,69 | 0,7 | 0,526 | | | |
| 100 g | 0,26 | 0,55 | 0,5 | 0,42 | 0,63 | 0,472 | | | |





AVK 1000.5.LLS comparator repeatability results

| Mari | MEASUREMENTS | | | | | | | |
|--------|--------------|-----------|-----------|-----------|-----------|---------|--|--|
| Mass | 1 SD [µg] | 2 SD [µg] | 3 SD [µg] | 4 SD [µg] | 5 SD [µg] | SD [µg] | | |
| 1000 g | 0,46 | 0,32 | 0,37 | 0,3 | 0,48 | 0,386 | | |
| 500 g | 0,27 | 0,44 | 0,4 | 0,31 | 0,45 | 0,374 | | |
| 200 g | 0,39 | 0,32 | 0,25 | 0,45 | 0,3 | 0,342 | | |
| 100 g | 0,36 | 0,21 | 0,17 | 0,4 | 0,36 | 0,3 | | |





Installations and implementations

Vacuum comparator with the LOAD LOCK system AVK 1000.5Y.LLS

Installation site: **GUM Polish NMI**Warsaw, Polska





Special application of the vacuum comparator

