





### Introduction

#### POLSKIE CENTRUM AKREDYTACJI POLISH CENTRE FOR ACCREDITATION



Sygnatariusz EA MLA

#### **CERTYFIKAT AKREDYTACJI**

LABORATORIUM WZORCUJĄCEGO ACCREDITATION CERTIFICATE OF CALIBRATION LABORATORY

Nr AP 069

Potwierdza się, że: / This is to confirm that:

#### RADWAG WAGI ELEKTRONICZNE WITOLD LEWANDOWSKI

ul. Toruńska 5, 26-600 Radom CENTRUM METROLOGII, BADAŃ I CERTYFIKACJI -LABORATORIUM POMIAROWE ul. Starowiejska 17A, 26-600 Radom

spełnia wymagania normy PN-EN ISO/IEC 17025:2018-02

Akredytowana działalność jest określona w Zakresie Akredytacji Nr AP 069 Accredited activity is defined in the Scope of Accreditation No AP 069

Akredytacja pozostaje w mocy pod warunkiem przestrzegania wymagań jednostki akredytującej określonych w kontrakcie Nr AP 069 This accreditation remains in force provided the Laboratory observes the requirements of Accreditation Body defined in the Contract No AP 069

> Akredytacji udzielono dnia 05.11.2004 r. ditation was granted on 05.11.200





POLSKIEGO CENTRUM AKREDYTACJI

Domilla LUCYNA OLBORSKA

Warszawa, dnia 4 stycznia 2021 roku

### Centre for Metrology, Testing and Certification – Measuring Laboratory of RADWAG

has the accreditation of the Polish Centre for Accreditation within the scope of pipettes calibration from 2009 as the first laboratory in Poland.

In 2023 it obtained the accreditation of the new ISO 8665-6:2022-11 norm within the scope of the single and multi-channel pipettes calibration.

#### ZAKRES AKREDYTACJI LABORATORIUM WZORCUJACEGO

SCOPE OF ACCREDITATION FOR CALIBRATION LABORATORY Nr/No AP 069

> wydany przez / issued by POLSKÍE CÉNTRUM AKREDYTACJI 01-382 Warszawa, ul. Szczotkarska 42

> > Wydanie/Issue 20 z/of 11.12.2023

Pipety tłokowe jednokanalowe	do 1 µl	0,050 µl	S	Procedura wewnetrzna			
Pipety tłokowe wielokanałowe	(1 ÷ 2) µl	0,050 µl		PW 05			
3 A	(2 ÷ 5) µl	0,050 µl		w oparciu o			
	(5 ÷ 10) µI	(5 ÷ 10) µl 0,050 µl					
	(10 ÷ 20) µl	0,07 µl		THE CHARLES THE STATE OF THE ST			
	(20 ÷ 50) µI	0,17 µl					
	(60 ÷ 100) µI	0,27 µI					
	(100 ÷ 300) µl	0,53 µI					
	(300 ÷ 500) µl	1,3 µl					
	(600 + 1000) µI	2,6 µl					
	(1000 ÷ 2500) µl	8,0 µl					
	(2500 ÷ 5000) µl	13 µl					
	(5000 + 10000) µl	20 µl					



### Introduction

# Major requirements relating to piston pipettes are regulated by the international norm ISO 8655

ISO 8655-1

Major requirements relating to piston pipettes are regulated by the international norm ISO 8655

- Part 1: Terminology, general requirements and recommendations for use

ISO 8655-2

Piston vessels for volume measurement

Part 2: Piston pipettes

ISO 8655-6

Piston vessels for volume measurement

Part 6: Gravimetric methods for defining of measurement errors



### **Minimal requirements** for scales

## Full compliance with ISO 8655-2/2022

Nominal volume of apparatus under test $(V)$	Resolution (d)	Repeatability (s)a	Expanded uncertainty in use $U(k=2)^{a, b}$			
	mg	mg	mg			
0,5 μl ≤ V < 20 μl	0,001 <sup>c</sup> 0,01 <sup>d</sup>	0,006 <sup>c, e</sup> 0,03 <sup>d</sup>	0,012 <sup>c, e</sup> 0,06 <sup>d</sup>			
20 μl ≤ V < 200 μl	0,01	0,025	0,05			
$200 \mu l \le V \le 10 ml$	0,1	0,2	0,4			
$10 \text{ ml} < V \le 1 000 \text{ ml}$	1	2	4			
$1000\mathrm{ml} < V \le 2000\mathrm{ml}$	10	10	40			

Single-channel balance.

Multi-channel balance, only valid for multi-channel pipettes. Multi-channel balances of 0,01 mg readability may be used to test multi-channel pipettes with nominal volumes below 20 µl only if the expanded uncertainty in use is less than one-fourth of the maximum permissible systematic error of the apparatus.



# Metrological requirements of piston pipettes

Volume of the liquid according to **ISO 8655** is related to mass and the Z rate contingent on liquid temperature and air pressure:

$$V = m \cdot Z$$

where: m - mass multiplied by

z - conversion rate [mg/ml]

$$[ml] = [mg \cdot \frac{1}{\frac{mg}{ml}} \cdot \frac{\frac{ml}{mg}}{\frac{ml}{ml}} = mg \cdot \frac{ml}{mg} = ml]$$



# Ambient conditions for conducting the measurements

Before starting the procedure, make sure that ambient temperature, the temperature of pipettes, tips and liquid

oscillate within the limits of  $(20 - 25)^{\circ}$ C and is maintained during measurements within the limits of  $\pm 0.5^{\circ}$ C **Relative humidity** should be  $(50 - 75)^{\circ}$ .



## **Dedicated solutions** of RADWAG company









MYA 5Y.P Microbalance for Pipette Calibration

XA 5Y.M.A.P Microbalance

**PSW Professional Single** Weighing Workstation

**AP-12.5Y** Automatic Device for Multichannel Pipette Calibration



## Minimum weigh requirements

## Full compliance with ISO 8655-2/2022

A multi-channel balance can be used to measure the test volume delivered from all channels in parallel, by aspirating and emptying each channel at the same time and analysing the results of each channel individually. The following variation may be applied to calibrate multi-channel pipettes with more than 12 channels or with 4,5 mm cone distance: use every second cone of the pipette to calibrate the pipette in multiple sets. All channels shall be tested.

If the pipette has more than one row of channels, it may be tested one row at a time.

Each channel can be measured individually, one after another, with a single-channel balance. For this purpose, test liquid shall be aspirated by all channels together and collected from one channel at a time. For the measurement of channel 1, for example, the volume of channel 1 is delivered into the weighing vessel, while the volumes from all other channels are discarded.



## **AP-12 Automatic Device for Multichannel Pipette Calibration**



**AP-12.1.5Y** [d]= 0,001 mg

**AP-12.5Y** [d] = 0.01 mg







Why?

- Speed, process automation, cost reduction
- Compliance with ISO 8655-2/2022



For whom?

- Accredited calibration laboratories
- Calibration centres
- Pipette users
- Pipette manufacturers




What do you gain?

- Compliance with ISO 8655-2/2022
- Time saving The procedure time is reduced from 6 hours to 1 hour
- Calibration of single and multi-channel pipettes



## **Advantages**





Single and Multichannel Pipette Calibration

- Calibration of up to 12-channel fixed-volume and variable-volume pipettes
- Automated calibration

#### **Calibration Weighing Pan**



As standard equipment, a weighing pan for checking and calibrating the balance is included.

#### **4 Times Faster**



The pump draws water from four measuring vessels at the same time. The procedure time is reduced from 6 hours to 1 hour.





## **PSW Professional Single Weighing Workstation**



## **PSW**

### **Professional Single Weighing Workstation**

To assure metrological traceability and meet all requirements arising from the measuring equipment supervision, and for the sake of easier and faster calibration of pipettes, RADWAG provides an effective and ergonomic automatic pipette calibration by the operator using the complex pipette calibration station.

#### Intended use:

Calibration of automatic piston pipettes:

- of a fixed volume
- of a variable volume



# PSW Professional Single Weighing Workstation



## **PSW**

**Professional Single Weighing Workstation** 

The working station is composed of three modules:

- 1. Measuring module
- 2. Environmental module
- 3. Computational module



# PSW Professional Single Weighing Workstation

#### **AMBIENT CONDITIONS**

The station is equipped with a thermohygrobarometer that allows monitoring ambient conditions on an ongoing basis.



- Water temperature measurement probe
- 2. Probe for measuring temperature, humidity and pressure



## **Other** balances









MYA 5Y.P Microbalance for Pipette Calibration

XA 5Y.M.A.P Microbalance

**XA 5Y.A Analytical Balances** 

**XA 5Y Analytical Balances** 



## **Other** balances



**XA 5Y.A Analytical Balances** 



**XA 5Y Analytical Balances** 

XA17 - Adapter for Pipette

Calibration



XA100 - Adapter for Pipettes calibration 100 ml XA 5Y.A



Nominal volume of apparatus under test V	Resolution (d) mg	Repeatability (s) <sup>3</sup> mg	Recommended balance										
			Microbalances				Analytical balances						
0,5 μl ≤ V < 20 μl	0,001 <sup>c</sup> 0,01 <sup>d</sup>	0,006 <sup>c, e</sup> 0,03 <sup>d</sup>		XA 6/21.5Y.M.A XA 6/21.5Y.M	XA 21.5Y.M.A XA 21.5Y.M	XA 21/52.5Y.M.A XA 21/52.5Y.M	XA 53.5Y.M.A XA 53.5Y.M						
20 µl ≤ V < 200 µl	0,01	0,025	MYA 21.5Y					XA 52.5Y	XA 110.5Y	XA 82/220.5Y.A XA 82/220.5Y		XA 120/250.5Y.A XA 120/250.5Y	
200 µl ≤ V ≤ 10 ml	0,1	0,2											
Adapter for pipette calibration		MY 11	XA 11	XA 11	XA 11	XA 11	XA 17	XA 17	XA 17	XA 100	XA 17	XA 100	
Maximum vessel capacity		11 ml	11 ml	11 ml	11 ml	11 ml	17 ml	17 ml	17 ml	100 ml	17 ml	100 ml	
Weighing vessel mass		~ 10 g	~ 10 g	~ 10 g	~ 10 g	~ 10 g	~ 11 g	~ 11 g	~ 11 g	~ 65 g	~ 11 g	~ 65 g	

