

PROFINET

Communication Protocol:
MW-01-A Mass Converter

SOFTWARE MANUAL

ITKP-40-01-12-21-EN



RADWAG
RADWAG WAGI ELEKTRONICZNE
ZAAWANSOWANE TECHNOLOGIE WAGOWE

DECEMBER 2021

CONTENTS

1. MASS CONVERTER SETTINGS AND CONFIGURATION	4
2. DATA STRUCTURE	4
2.1. Input Address	4
2.2.1 Input Registers	4
2.2. Output Address	7
2.2.1 Output Registers	7
3. CONFIGURATION OF THE PROFINET MODULE IN TIA PORTAL V16 ENVIRONMENT	10
3.1. Import GSD	10
3.2. Module Configuration	12
4. DIAGNOSTICS APP	16

1. MASS CONVERTER SETTINGS AND CONFIGURATION

To configure MW-01-A mass converter settings for communication via PROFINET protocol, run **MwManager** PC software and go to **<Parameters / Communication / Additional modules>**. For detailed description of configuration read **MWManager** user manual.

2. DATA STRUCTURE

All registers are 2-byte type (WORD). Floating point data (such as mass and tare) are stored in two consecutive registers and are FLOAT-type data. If the first register consists of two AB bytes and the other of two CD bytes, then FLOAT value is HEX ABCD,

2.1. Input Address

Input variables:

Variable	Offset	Length [WORD]	Data type
Platform mass	0	2	float
Platform tare	4	2	float
Platform unit	8	1	word
Platform status	10	1	word
Platform LO threshold	12	2	float
Process status	64	1	word
Input status	66	1	word
Min	68	2	float
Max	72	2	float
Fast dosing threshold	76	2	float
Slow dosing threshold	80	2	float
Adjustment status	100	1	word

2.2.1 Input Registers

Platform mass – returns platform net mass in a current unit.

Example:

The read register with offset 0 has a hex value of 0x43E28000, after converting to float we get 453.0 which is the current indication of the load mass.

Platform tare – returns platform tare in an adjustment unit.

Platform unit – determines a current (displayed) net mass unit of the platform.

Unit bits	
0	Gram [g]
1	Kilogram [kg]
2	Carat [ct]
3	Pound [lb]
4	Ounce [oz]
5	Newton [N]

Example:

Read HEX value: 0x02. Binary form:

B1/7	B1/6	B1/5	B1/4	B1/3	B1/2	B1/1	B1/0	B0/7	B0/6	B0/5	B0/4	B0/3	B0/2	B0/1	B0/0
0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0

The mass unit is kilogram [kg].

Platform status – determines status of a weighing platform.

Status bits	
0	Measurement correct (the scale does not report any error)
1	Measurement stable
2	Scale indicates zero
3	Scale tared
4	Scale in II weighing range
5	Scale in III weighing range
6	Scale reports NULL error
7	Scale reports LH error
8	Scale reports FULL error

Example:

Read HEX value: 0x13

B1/7	B1/6	B1/5	B1/4	B1/3	B1/2	B1/1	B1/0	B0/7	B0/6	B0/5	B0/4	B0/3	B0/2	B0/1	B0/0
0	0	0	0	0	0	0	0	0	0	0	1	0	0	1	1

The scale does not report any error, measurement stable in weighing range II.

LO threshold – returns value of platform's **LO** threshold in an adjustment unit.

Process status – determines status of the dosing/formulation process:

HEX value	Description
0x00	Process disabled
0x01	Process started
0x02	Process aborted
0x03	Process finished

Input state – bitmask of mass converter inputs. Three least significant bits represent the mass converter inputs state.

Example:

Read HEX value: 0x0005

B1/7	B1/6	B1/5	B1/4	B1/3	B1/2	B1/1	B1/0	B0/7	B0/6	B0/5	B0/4	B0/3	B0/2	B0/1	B0/0
0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	1

Inputs number 1 and 3 of the weighing terminal take HI state.

MIN – returns **MIN** threshold value in an adjustment unit.

MAX – returns **MAX** threshold value in an adjustment unit.

Fast dosing threshold - returns fast dosing threshold value in an adjustment unit.

Slow dosing threshold - returns slow dosing threshold value in an adjustment unit.

Adjustment status – determines adjustment process status:

HEX value	Description
0x00	Process finished correctly
0x01	Start mass/adjustment coefficient determination in progress
0x02	Range exceeded
0x03	Time exceeded
0x04	Process aborted
0x05	Awaiting for data

2.2. Output Address

Output variables list:

Variable	Offset	Length [WORD]	Data type
Basic command	0	1	word
Command with parameter	2	1	word
Tare	6	2	float
LO threshold	10	2	float
Outputs state	14	1	word
Min	16	2	float
Max	20	2	float
Fast dosing threshold	24	2	float
Slow dosing threshold	28	2	float

2.2.1 Output Registers

Basic command – writing a respective value triggers direct performance of a given task, see the table:


Bit No.	Operation
0	Zero platform
1	Tare platform
5	Process start
6	Process stop

Example:

Writing the register with value 0x02



B1/7	B1/6	B1/5	B1/4	B1/3	B1/2	B1/1	B1/0	B0/7	B0/6	B0/5	B0/4	B0/3	B0/2	B0/1	B0/0
0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0

This causes scale taring.

	<p><i>A command is executed once upon detecting that its bit has been set. If the command is to be executed more than once, it is necessary to zero the bit first, and reset it to the required value next.</i></p>
---	--

Complex command (with parameter) – setting a respective value triggers performance of a given task, see the table:

Bit No.	Operation
0	Setting tare value for a platform
1	Setting LO threshold value for a platform
2	Setting outputs state
3	Setting MIN threshold value
4	Setting MAX threshold value
5	Setting fast dosing threshold
6	Setting slow dosing threshold

	<i>Complex command requires setting a respective parameter (offset from 6 to 28 – refer to output registers table)</i>
	<i>A command with a parameter is executed once when its bit setting is detected. If the command is to be executed more than once, it is necessary to zero the bit first, and reset it to the required value next.</i>

Example:

Sending tare of 1.0 value to the scale.

Performance of the command requires writing of 2 registers:

offset 2 – command with a parameter - value 0x01 – i.e. tare setting,

offset 6 – tare value in float format - 1.0.

Tare – complex command parameter: tare value (in an adjustment unit).

LO threshold – complex command parameter: LO threshold value (in an adjustment unit).

Output state – complex command parameter: state of mass converter outputs

Example:

Setting high state for output 1 and 3 of the mass converter.

Output mask:

B1/7	B1/6	B1/5	B1/4	B1/3	B1/2	B1/1	B1/0	B0/7	B0/6	B0/5	B0/4	B0/3	B0/2	B0/1	B0/0
0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	1

Upon conversion to HEX the result is 0x05

Performance of the command requires record of 2 registers:

offset 2 – command with parameter - value 0x04 – i.e. output state record,

offset 14 – output mask 0x05.

As a result, outputs number 1 and 3 take high state.

MIN – complex command parameter: MIN threshold value (in an adjustment unit).


MAX – complex command parameter: MAX threshold value (in an adjustment unit).

Fast dosing threshold - complex command parameter: fast dosing threshold value (in adjustment unit).

Slow dosing threshold - complex command parameter: fine dosing threshold value (in adjustment unit).

3. CONFIGURATION OF THE PROFINET MODULE IN TIA PORTAL V16 ENVIRONMENT


Start operation in the environment by setting up a new project, where topology of PROFINET network with MASTER controller will be specified, in this case the MASTER controller is SIEMENS controller of S7-1200 series.




3.1. Import GSD

Using a delivered configuration file (GSD), add a new device to the environment.


Open OPTIONS tab, next click MANAGE GENERAL STATION DESCRIPTION FILES (GSD) entry and select access path to the GSD file.



With the file successfully added, the ABIC-PRT module is displayed on the list of devices.



Now, it is possible to make a network comprising a MASTER controller and a newly added SLAVE module.



3.2. Module Configuration

At this stage, it is necessary to create a network consisting of MASTER controller and SLAVE device (MW-01-A mass converter). With the power supply on, you can use ACCESSIBLE DEVICES function to search for a given device in the environment. The list should display both the MASTER and the SLAVE.

Type of the PG/PC interface:

PG/PC interface:

Accessible nodes of the selected interface:

Device	Device type	Interface type	Address	MAC address
pawelk	SIMATIC-PC	PN/IE	10.10.3.145	4C-ED-FB-44-C0-31
mw01-profinet	RT Migration (FW 1.xx)	PN/IE	10.10.8.100	00-30-11-34-44-E6
plc_1	CPU 1211C AC/DC/Rly	PN/IE	10.10.8.222	E0-DC-A0-CF-59-E0
tk-ka-cperczyk-m	SIMATIC-PC	PN/IE	10.10.22.21	F8-32-E4-A0-BF-29

Flash LED


Start search

Online status information: Display only error messages


- Found accessible device mw01-profinet
- Scan completed. 4 devices found.
- Retrieving device information...
- Scan and information retrieval completed.



Show Cancel

Next, specify the IP address of the module and its name in PROFINET network. To do it, select the module, in PROPERTIES tab find PROFINET INTERFACE submenu, click it and enter the IP address and name. Make sure that these settings are accordant with parameters set in the MW-01 mass converter. Remember that the SLAVE IP address and MASTER address must come from the same subnet.




Next, configure the module. First, specify input and output registers size, define their start addresses. From the list of available INPUT and OUTPUT modules select such modules as presented in the picture below. Maximum size of input data is 102 bytes, maximum size of output data is 52 bytes. The project uses the default starting addresses: 256 for INPUT and 256 for OUTPUT with a size of 102 bytes for input data and 32 bytes for output data, as presented below.






At this stage, it is possible to upload the hardware configuration into the controller.




Upon successful compilation and code reading, the MASTER and SLAVE shall establish communication. It can be verified via ONLINE connection. The result should be as presented below:



Now, proceed to the process of program code making.

4. DIAGNOSTICS APP

Start creating the app by defining names of symbolic input and output registers. For this purpose, use a branch of the PLC TAGS tree. For the purpose of this example, tag tables have been created as presented below:



PROFINET input and output registers are specified in TAG_ProfiNET_IN and TAG_ProfiNET_OUT tables. The below screenshots present symbolic names and addresses:

The screenshot shows the 'TAG_ProfiNET_IN' tag table in a software interface. The table has columns for Name, Data type, Address, Retain, Acces..., Writa..., Visibi..., and Comment. The table contains 12 rows of data, each with a small icon in the first column.

	Name	Data type	Address	Retain	Acces...	Writa...	Visibi...	Comment
1	AnyBus_Platform_mass	DWord	%ID256	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
2	AnyBus_Platform_tare	DWord	%ID260	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
3	AnyBus_Platform_unit	Word	%IW264	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
4	AnyBus_Platform_status	Word	%IW266	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
5	AnyBus_Platform_LO_threshold	DWord	%ID268	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
6	AnyBus_Process_status	Word	%IW320	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
7	AnyBus_Inputs_status	Word	%IW322	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
8	AnyBus_Min	DWord	%ID324	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9	AnyBus_Max	DWord	%ID328	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10	AnyBus_Fast_dosing_threshold	DWord	%ID332	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
11	AnyBus_Slow_dosing_threshold	DWord	%ID336	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
12	AnyBus_Adjustment_status	Word	%IW356	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	


Project_08_Profinet_MW_01_New_instrukcja ▶ PLC_1 [CPU 1211C AC/DC/Rly] ▶ PLC tags ▶ TAG_ProfiNET_OUT [9]

Tags User constants

TAG_ProfiNET_OUT

	Name	Data type	Address	Retain	Acces...	Writa...	Visibl...	Comment
1	AnyBus_Command	Word	%QW256	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
2	AnyBus_Command_with_para...	Word	%QW258	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
3	AnyBus_Tare	DWord	%QD262	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
4	AnyBus_LO_threshold	DWord	%QD266	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
5	AnyBus_Output_state	Word	%QW270	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
6	AnyBus_Set_Min	DWord	%QD272	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
7	AnyBus_Set_Max	DWord	%QD276	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
8	AnyBus_Set_Fast_dosing_thres...	DWord	%QD280	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9	AnyBus_Set_Slow_dosing_thre...	DWord	%QD284	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10	<Add new>			<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	

So as not to work directly on the physical inputs/outputs of the module, it is recommended to create data blocks comprising representations of these registers, and “copy” values between them. For this purpose, define data blocks as follows:




ProfiNET_Input and ProfiNET_Output blocks represent input/output registers of the MW-01-A mass converter's PROFINET module. See the screenshots below:

Name	Data type	Offset	Start value	Retain	Accessible f...	Writa...	Visible in ...	Setpoint	Comment
1	Static								
2	mass	Real	0.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3	tare	Real	4.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4	unit	Word	8.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5	platform_status	Word	10.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6	LO	Real	12.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7	process_status	Word	16.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8	inputs	Word	18.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
9	min	Real	20.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
10	max	Real	24.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
11	fast_dosing_threshold	Real	28.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
12	slow_dosing_threshold	Real	32.0	0.0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
13	adjustment_status	Word	36.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	

Name	Data type	Offset	Start value	Retain	Accessible f...	Writa...	Visible in ...	Setpoint	Comment
1	Static								
2	Command	Word	0.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3	Command_with_parameter	Word	2.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4	Platform	Word	4.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5	Tare	DWord	6.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6	LO_threshold	DWord	10.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7	Output_state	Word	14.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
8	Min	DWord	16.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
9	Max	DWord	20.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10	Fast_dosing_threshold	DWord	24.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
11	Slow_dosing_threshold	DWord	28.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
12	Adjustment_weight_mass	DWord	32.0	16#0	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	

A MOVE instruction, for example, can be used to rewrite values between the physical inputs/outputs of the module and the registers in the data blocks:



Upon compilation and upload of the program to the controller in data block, it is possible to read input registers (MONITOR ALL) and to record output registers (e.g. by changing START VALUE and LOAD START VALUES AS ACTUAL) of the SLAVE module.



RADWAG WAGI ELEKTRONICZNE
ZAAWANSOWANE TECHNOLOGIE WAGOWE

