

**FULL-FEATURED SOLUTION** 





# Pipettes Calibration Workstation



### **FULL-FEATURED** Workstation

Anti-vibration weighing table

Sink and container for used liquid

Temperature, humidity and pressure probes

> Distilled water temperature probe

PC with keyboard, mouse and LCD

Pipettes PC Software:

Calbration of adjustable volume pipettes(channels quantity:1, 8, 12)

Automated measurement procedures

Pipettes database and Calibration Results database

Result calculation based on: mean pipette volume (channel), systematic error e (error of accuracy A) random error CV (error of repeatability P)

## **Full-featured Solution for Pipettes Calibration Excellent ergonomics and precise measurement**

#### The Highest Measurement Accuracy

An integral part of the pipettes calibration workstation is either MYA 21.5Y.P microbalance, XA 52.5Y.M.A.P microbalance or XA 82/220.5Y.A analytical balance, both featuring a special adapter for pipettes calibration. The devices are characteristic for their evaporation rings limiting pipette calibration errors to the minimum. After dimantling od the adapter, each device can be used for standard weighing processes.

#### **High Resistance to Ground Vibrations**

An in-built anti-vibration table featuring granite top prevents potential ground vibrations from being transferred to the microbalance. With this solution stabilization time shortens considerably. Separation of chipboard-made table components from anti-vibration frame blocks vibrations being an effect of user operation.

#### **Ambient Conditions Monitoring**

Set of probes, used to measure humidity, atmospheric pressure and temperature of air and distilled water, uninterruptedly monitors ambient conditions. This guarantees absolute reliability of carried out measurements

#### **Computer Support for Pipettes Calibration Process**

Pipettes PC Software facilitates automation of pipettes calibration process carried out using gravimetric method, accordant with an international standard, ISO 8655-6. In compliance with ISO 10012 Standard, the software enables generating reports on calibration, archiving measurements' results and overall pipettes' management.



MYA 21.5Y.P Microbalance

Max capacity [Max]: 21 g

Radability [d]:

Minimum Sample Weight by USP: 2 mg

Weighing pan ø 26 mm + adapter for pipettes calibration

> Automatic door opening system

colour touchscreen

2-point internal adjustment

Interfaces: 2×USB-A, USB-C, HDMI, Wi-Fi®, Ethernet, Hotspot



XA 52.5Y.M.A.P Microbalance

Max capacity [Max]: 52 g

> Radability [d]: 5 μg

Minimum Sample Weight by USP: 4.4 mg

Weighing pan ø 26 mm + automatic adapter for pipettes calibration

> Automatic door opening system

colour touchscreen

2-point internal adjustment

Interfaces: 2×USB-A, USB-C, HDMI, Wi-Fi®, Ethernet, Hotspot



XA 82/220.5Y.A **Analytical Balance** 

Max capacity [Max]: 82 g / 220 g

Radability [d]: 0.01 mg / 0.1 mg

Minimum Sample Weight by USP: 10 mg

Weighing pan ø 90 mm + adapter for pipettes calibration

> Automatic door opening system

colour touchscreen

2-point internal adjustment

Interfaces: 2×USB-A, USB-C, HDMI, Wi-Fi®, Ethernet, Hotspot

Wi-Fi® is a registered trademark of Wi-Fi Alliance.



**SINGLE:** full-featured workstation + a selected balance (MYA 21.5Y.P or XA 52.5Y.M.A.P or XA 82/220.5Y.A)

COMBO: full-featured dual workstation + MYA 21.5Y.P microbalance + a selected balance (XA 52.5Y.M.A.P or XA 82/220.5Y.A)