

User Manual

WWG-2 series large-size display

Manual number:
ITKU-07-03-08-14-EN



- **WWG-2/3**
- **WWG-2/4**



BALANCES AND SCALES

RADWAG, 26 – 600 Radom, Bracka 28
Phone +48 48 38 48 800, phone/fax: +48 48 385 00 10
Sales department +48 48 366 80 06
www.radwag.com

AUGUST 2014

Table of contents

| | | |
|----|-------------------------------------|---|
| 1. | INTENDED USE | 4 |
| 2. | TECHNICAL SPECIFICATION | 4 |
| 3. | DIMENSIONS..... | 5 |
| 4. | MOUNTING..... | 5 |
| 5. | START-UP | 5 |
| 6. | SETTING PARAMETERS..... | 6 |
| 7. | COOPERATION WITH AN INDICATOR | 6 |
| 8. | COMMUNICATION PROTOCOL..... | 8 |

1. INTENDED USE

WWG-2 large-size display comprises six backlit LCDs with digits size of 4,5" (11,2 cm). The display is housed in a stainless steel casing. It is equipped with a mounting bracket allowing either for placing the display on a table or for hanging it on a wall.


Communication with the indicator is established by means of **RS232** interface. **WWG-2/3** design is intended for operation with **PUE C/31H** indicator. **WWG-2/4** design is intended for operation with both, **PUE C/41H** indicator and **PUE 5** terminal. The display uses white diode backlight technology. It is powered with **230V AC**.

WWG displays mass measured by means of industrial scales. It is especially useful when measuring is performed in severe ambient conditions, i.e. when bright, intense, distant light (including sun) disrupts clear visibility.


2. TECHNICAL SPECIFICATION

| | WWG-2/3 | WWG-2/4 |
|------------------------------------|-------------------------|-----------------|
| Housing | Stainless steel | |
| Digits size | 4,5" (11,2cm) | |
| Intended for indicator | PUE C/31H | PUE C41H, PUE 5 |
| Communication | RS232 | |
| Means of connecting | interface | |
| Signal cable length | 3m (optionally max 10m) | |
| Power supply | 230VAC | |
| IP rate | IP66/67 | |
| Working temperature | od -20°C do +40°C | |
| Dimensions (mounting kit excluded) | 635x195x48 | |
| Mass (nett value) | 5,5kg | |

3. DIMENSIONS



4. MOUNTING



5. START-UP

Depending on cable type, the display shall be connected to RS232 interface of PUE C/31H and PUE C/41H indicator or PUE 5 terminal. Next it must be connected to mains. Upon establishing communication the display is backlit, the mass displayed.

Caution:

For correct cooperation of an indicator and a display it is necessary to set the display parameters in accordance with point 6.

6. SETTING PARAMETERS

PUE C/31H parameters:

| | Parameter name and number | Parameter value |
|-------------------------|---------------------------|-----------------|
| Printout type | 2.1. Pr_n | CntA |
| Baud rate | 2.3. bAud | 2400 - 38400 |
| Transmission parameters | 2.4. S_rS | 8d1SnP |

PUE C41H parameters:

| | Parameter name and number | Parameter value |
|-------------------------|---------------------------|--------------------|
| Send Data Port | P 3.5.1 DISP. PORT | RS232(1), RS232(2) |
| Additional display type | P 3.5.2 DISP. TYPE | WWG |
| Baud rate | P 2.2-3.1 BAUD RATE | 2400 - 38400 |
| Transmission parameters | P 2.2.2 DATA BITS | 8 |
| | P 2.2.3 PARITY BITS | NONE |
| | P 2.2.4 STOP BITS | 1 |

7. COOPERATION WITH AN INDICATOR

RS232 pins marking for WWG-2 large-size display boards:




Diagram of cable connecting PUE C/31H with RS232 port of WWG-2/3 large-size display:





Diagram of cable connecting PUE C/41H, PUE 5 with RS232 port of WWG-2/4 large-size display:



8. COMMUNICATION PROTOCOL

The display supports frame formats sent in a form of continuous transmission.

Frame format:

| | | | | | | | | | | | | |
|---|---|-------|------------------|-------|-----------|------|-------|------|----|----|----|----|
| 1 | 2 | 3 | 4 | 5 | 6 | 7-15 | 16 | 17 | 18 | 19 | 20 | 21 |
| S | I | space | Stability marker | space | character | mass | space | unit | | | CR | LF |

Stability marker

[space] if stable

[?] if unstable

[^] if error of exceeding range to + occurs

[v] if error of exceeding range to - occurs

Character

[space] for positive values

[-] for negative values

Mass

9 characters with decimal point and right justification

Unit

3 characters with left justification

Command

3 characters with left justification

**MANUFACTURER
OF ELECTRONIC WEIGHING
INSTRUMENTS**



RADWAG BALANCES AND SCALES
POLAND, 26-600 Radom, Bracka 28

Phone. +48 48 38 48 800, fax. + 48 48 385 00 10

Sales department + 48 48 366 80 06

export@radwag.com

www.radwag.com

