

BATH GEL – SHOWER GEL dry matter content determination

High-quality cosmetic products for personal hygiene, such as gels, must fulfill their basic roles that is washing, moisturizing, regenerating and relaxing. The efficiency of each of these stages is the result of selection of suitable ingredients of a proper mass and concentration in the gel. Real amounts of each ingredient added to the mixture must be supervised, which eventually determines the quality of the gel and product acceptance by the user. One of engineering parameters of numerous cosmetic products is dry mass, that is part of the sample that remains after removing all volatile ingredients. The method of measuring the dry matter to be used in testing must guarantee accuracy and highly precise results, which can be achieved with the use of MA/R and MA/X2 moisture analyzer by Radwag.



The application note includes basic information for validation of the bath gel – shower gel drying method with the use of MA/R and MA/X2 moisture analyzers series by Radwag Wagi Elektroniczne. The application note may be the basis for elaborating own drying method with special regard to distinctive features of the product in question.



Bath gel, shower gel – dry matter content determination

The method with the use of IR radiation Metrology, Research and Certification Center, Radwag Wagi Elektroniczne, Poland Toruńska 5, 26-600 Radom, Poland +48 48 386 60 00, e-mail: office@radwag.com, <u>www.radwag.com</u>

TERMS

ACCURACY of determining water / dry matter content is the difference between the result of the water / dry matter content received in the moisture analyzer method and the result of the water / dry matter content received while drying the same sample through a reference method.

PRECISION is a degree of compliance between independent results of the test, received in specific conditions. The measure of precision is a standard deviation from a series of several measurements.

REFERENCE METHOD

The reference method parameters are usually specified in standards or other discipline-specific documents as the socalled guides. If such documents are unavailable, the drying temperature that does not cause the sample to change colors is used. Such an approach applies to previously dehydrated products and raw products.

SAMPLE PREPARATION

Before testing, sample must be stored in sealed containers. Fluid and semi-fluid products must be mixed before testing.

ACCESSORIES

Dryer, quartz sand, weighing vessels with a lid, rods, AS 220.X2 balance, laboratory spoon.

METHOD DESCRIPTION

Weigh glass vessels with a glass rod and pre-dried quartz sand in the amount of ca. 15 g. Mix the sample with a mass of ca. 5 g and put into the glass weighing vessels on pre-dried quartz sand. Mix the sample with sand by means of the glass rod that must be left in the vessel. Use of sand as a foundation is aimed at elimination creation of the shell on the surface of the sample in question. Weigh vessels again and specify the real mass of the sampl in question with the use of the balance whose weighing accuracy is 0.1 mg (AS 220.X2). Place weighing vessels with the sample and lids in the temperature-controlled laboratory dryer. Dry samples at the temperature of 105°C for 3 hours. After this period, remove vessels and put into the desiccator until they cool down and weigh afterwards. Place samples in the laboratory dryer again and keep on drying them for 30 minutes. Cool them down and weigh again. Repeat the procedure until you obtain a stable sample mass or record the sample mass growth after drying. Calculate the dry matter content as a quotient of the post-drying sample mass (m₂) and pre-drying sample mass (m₁)

RESULTS

Sample name	BATH GEL	SHOWER GEL
Dry matter content (%)	16.43	16.09
Standard deviation (%)	0.12	0.08

BATH GEL AND SHOWER GEL – DRY MATTER CONTENT DETERMINATION WITH THE MOISTURE ANALYZER

The water content testing with the use of the moisture analyzer (IR radiation) entails two phenomena: convection and radiation. The sample temperature rises from outer layers to the bottom of the sample. The temperature gradient in the sample structure minimizes through optimization of the thickness of the dried sample and drying temperature.

SAMPLE PREPARATION

Before testing, samples must be stored in sealed packaging. Fluid and semi-fluid products must be mixed before testing.

ACCESSORIES

MA/R or MA/X2 moisture analyzer, glass weighing vessels with a lid, quartz sand, laboratory spoon, pipette.

METHOD DESCRIPTION

Set drying parameters presented below. Dose the sample with a mass of ca. $1.5 \div 2$ g by means of the pipette onto pre-dried quartz sand. Lock the drying chamber manually or automatically. Note: it is possible to dose the sample straight into the moisture analyzer weighing pan, but the analysis time would be around twice longer.

DRYING PARAMETERS / RESULTS

Sample name	BATH GEL	SHOWER GEL
Drying profile	Standard	
Drying temperature	130°C	120°C
Sample mass (g)	~ 1.5 ÷ 2	
End of analysis	Auto 1	
Dry matter content (%)	16.93	16.53
Standard deviation (%)	0.19	0.07
Analysis time \acute{x} (min)	~ 7	~ 8

ACCURACY OF THE MA/R ÷ MA/X2 METHOD

Sample name	BATH GEL	SHOWER GEL
Dry matter content (%) – Ref.	16.43 ± 0.12	16.09 ± 0.08
Dry matter content (%) – MA R/X2	16.93 ± 0.19	16.53 ± 0.10
Analysis accuracy (%)	0.50	0.44

RESERVATION

The method in question has been verified by the Research Laboratory, yet the results do not include factors arising from diversity of tested samples, operators' personal skills as well as measuring capability used by moisture analyzer users. For this reason Radwag shall not be held responsible for drying parameters but they can be used to elaborate own drying method.

