



More information on the website
radwag.com/it/info,w1,12Q

Suspended Self-Centring Weighing Pan for WAY 500.5Y.KO, WAY 1.5Y.KO, WAY 2.5Y.KO Manual Mass Comparator



The drawings, photos and graphics used are for illustrative purposes only.

Caratteristiche generali

The suspended self-centring weighing pan is designed specifically for the WAY 5Y.KO manual mass comparators.

Installed It is installed during the device's production stage in place of the standard weighing pan, so it should be ordered when ordering a comparator (it is not sold as a separate product). The weighing pan has the following functions: placing several weights, which are the resultant of the reference standard, off-centre, positioning the centre of gravity of the test and reference standards always in the same place.



Compatibile con

WAY 5Y.K0 Comparatore di massa manuale