



More information on the website
radwag.com/en/info,w1,6NX

Cable PT0453

WX-013-0238



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

Compatible with (Additional Fee)

XA 5Y.M.A Microbalance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance
Automatic feeders
XA 5Y.M Microbalance
Receipt Printer