



More information on the website  
[radwag.com/en/info,w1,A7T](http://radwag.com/en/info,w1,A7T)

# 2 $\mu$ l ÷ 20 $\mu$ l - RW8-102-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
Volume	2 $\mu$ l ÷ 20 $\mu$ l
Systematic error	3.0% / 0.9%
Random error	2.0% / 0.4%
Tips	300 $\mu$ l



## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
 Microbalance for Pipette Calibration MYA 5Y.P

Automatic Balances for. Pipettes Calibration  
 Pipette Tips