

AS 3100.X7 Analytical Balance





The drawings, photos and graphics used are for illustrative purposes only.

Functions

Q	Autotest		Dosing	- <u>0K</u> +	Plus/Minus Control	%	Percent Weighing
***	Parts counting	MAX	Peak hold		Formulation	7	Newton unit measurement
<u>.al</u>	Statistics	- <u>0K</u> +	Checkweighing	4	IR sensors	8	Under-pan weighing
GLP	GLP Procedures		Animal weighing	ρ	Density determination		Ambient conditions monitoring
4	Replaceable unit	SQC	Statistical Quality Control		ALIBI Memory	141	Drying modes
2))))	Samples drying	♦	Moisture content analysis	- ` ↓- %D	Dry mass determination	Ш	Mass for titrator
	W. F.						

Datasheet

Metrological parameters	
Maximum capacity [Max]	3100 g
Minimum load	- mg
Readability [d]	1 mg
Verification unit [e]	-
Tare range	-3.1 kg
Standard repeatability [5% Max]	0.5 mg
Standard repeatability [Max]	0.6 mg
Standard minimum weight (USP)	1 g
Standard minimum weight (U=1%, k=2)	100 mg
Permissible repeatability [5% Max]	0.8 mg
Permissible repeatability [Max]	1 mg
Linearity	±4 mg
Stabilization time	2 s
Adjustment	internal (automatic)
OIML Class	-
Dhysical nevernators	
Physical parameters	comi sustamentia. LauraICENICINIC
Leveling system	semi-automatic – LevelSENSING
Display	7" graphic colour touchscreen
Weighing chamber doors	manual Balance, weighing pan, weighing pan shield, centring ring, bottom
Delivery components	cover, power supply.
Weighing chamber dimensions	190×190×222 mm
Weighing pan dimensions	ø90 mm (open-work pan)
Packaging dimensions	490×400×520 mm
Net weight	7.3 kg
Gross weight	9.3 kg
Construction	
Protection class	IP 43
Commonants and aufturan	
Components and software Database capacity	7
Database Capacity	,
Features of use	
Touch-free operation	2 IR Sensors
Communication interface	
Communication interface	RS232¹, 2×USB-A (interchangeable), USB-B, Wi-Fi, Ethernet
Electrical parameters	Adapter 100 040V 40 50/60U- 0 64 40V 50 4 04
Power supply	Adapter: 100 – 240V AC 50/60Hz 0.6A; 12V DC 1.2A Balance: 12 – 15V DC 0.8A max
Power consumption max.	4 W
Environmental conditions	
Operating temperature	+10 - +40 °C
Ambient conditions monitoring (option)	THBR 2.0 System, THBR BOX, THB P, THB W, THB S
, and the conditions monitoring (option)	THER 2.0 System, THER BOX, THET, THE W, THE S

Environmental conditions

Relative humidity

40% - 80%

Repeatability is expressed as a standard deviation from 10 weighing cycles.

Stabilization time depends on the ambient conditions and the dynamics of weighing pan loading; specified for FAST profile.

¹ Barcode scanners, available as weighing instrument accessory, communicate with the instrument via RS232 interface exclusively.

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.



Accessories

Antivibration Tables
Holders for laboratory flasks
Power Adapters
Cigarette lighter receptacle power supply cables
Density determination KIT
USB cable (scale - printer)
Professional Weighing Tables
Barcode scanners
Holders for test tubes and filters
Workstation for Pipettes Calibration
RS 232, RS 485 cables

THBR 2.0 System - Ambient Conditions Monitoring Displays
Protective cover for balances
Weighing dishes
Antistatic ionizer
Receipt Printer
RS 232, RS 485 cables
Additional modules
Under-pan weighing
RS 232 cables (scale - printer)
RS 232 – RS 485 Converter

Software

- RAD Key [WX-010-0005]
- R-Lab [WX-010-0080]
- RADWAG Development Studio [WX-010-0104]
- Alibi Reader [WX-010-0114]
- Scale Editor 2.1 [WX-010-0173]

Device dimensions

