



More information on the website  
[radwag.com/en/info,w1,C80](http://radwag.com/en/info,w1,C80)

# 10 $\mu$ l ÷ 100 $\mu$ l - RW8-104-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
Volume	10 $\mu$ l ÷ 100 $\mu$ l
Systematic error	3.0% / 0.8%
Random error	1.5% / 0.15%
Tips	300 $\mu$ l



## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
Microbalance for Pipette Calibration MYA 5Y.P

Automatic Balances for. Pipettes Calibration  
Pipette Tips