



More information on the website  
[radwag.com/us/info,w1,C80](https://radwag.com/us/info,w1,C80)

# 10 $\mu$ l ÷ 100 $\mu$ l - RW8-104-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Volume	10 – 100 $\mu$ l
Systematic error	3,0% / 0,8%
Random error	1,5% / 0,15%
Tips	300 $\mu$ l



## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
MYA 5Y.P Microbalance for pipettes calibration

Automatic Balances for. Pipettes Calibration  
Pipette Tips