



More information on the website
radwag.com/us/info,w1,0C3

PT0455 Cable



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0455 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a RADWAG RTP printer. Its standard length is 1,8m.

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration
UMA 5Y.F Automatic Weighing System
5Y.PM Precision balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
Susceptometer - Magnetization
Automatic Balances for. Pipettes Calibration
WAY 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
UYA 5Y Ultra-Microbalance
AK-4 Automatic Mass Comparator
AKM-2 Automatic Mass Comparator
Analytical Balance for Filters XA 5Y.F
UMA Automatic Mass Comparator

Density measurement
MYA 5Y Microbalance
5Y.PM.KB Manual Mass Comparator
HRP 5Y.KO Mass Comparator
Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
HRP 5Y.KB Mass Comparator
RMC F Robotic Weighing System
RB F Robotic Weighing System
XA 5Y.M Microbalance
UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance for stents weighing

