



More information on the website  
[radwag.com/us/info,w1,QLH](http://radwag.com/us/info,w1,QLH)

# 5 µl ÷ 50 µl - RW8-103-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Volume	5 – 50 µl
Systematic error	2,0% / 0,6%
Random error	2,0% / 0,3%
Tips	300 µl



## Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration

Automatic Balances for. Pipettes Calibration  
 Pipette Tips