



More information on the website  
[radwag.com/tr/info,w1,QNI](http://radwag.com/tr/info,w1,QNI)

## Voluminous sample container



The drawings, photos and graphics used are for illustrative purposes only.

### Genel Açıklama

The set is an ideal solution for drying large volume samples, e.g. textiles, or products that placed on the weighing pan tend to increase their volume. The container allows to enclose the sample within a limited space wherein the possibility to determine the moisture content value remains unaffected.

In the case of large-volume samples, the kit prevents the sample from coming into contact with solid elements of the moisture analyzer (such as the heating element or its cover) and protects it from damage. Closing a large-volume, low-density sample in the container also reduces the risk of errors resulting from the effect of the sample being raised by warm air.

The accessory gives the user full control of the sample and an error-free result.

The set consists of 2 components made of steel. It is compatible with all moisture analyzers.

### İle uyumlu

MA X7.IC.A Nem Tayin Cihazı  
MA 5Y Nem Tayin Cihazı  
MA X2.IC.A Nem Tayin Cihazı

MA X7.A Nem Tayin Cihazı  
MA 5Y.IC Nem Tayin Cihazı  
MA R Nem Tayin Cihazı

