




More information on the website
radwag.com/us/info,w1,Q01

Fingerprint Reader



The drawings, photos and graphics used are for illustrative purposes only.

Overview

An accessory for the ELLIPSIS laboratory balances, designed to read fingerprints in order to log into the balance. Simply touch the reader to log in. The fingerprint reader is connected to the ELLIPSIS balance via the USB port.



Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration
AK-5/100.STENT Automatic Balance for Stents
UMA 5Y.F Automatic Weighing System
MA 5Y Moisture Analyzer
5Y.PM Precision balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
PUE CY10 Weighing Terminal

UMA Automatic Mass Comparator
APP 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
Density measurement
MYA 5Y Microbalance
XA 5Y Analytical Balance
5Y.PM.KB Manual Mass Comparator
HRP 5Y.KO Mass Comparator
PS 5Y Precision Balance

CY10 Multifunctional Scale
Susceptometer - Magnetization
Automatic Balances for. Pipettes Calibration
XA 5Y.A Analytical Balance
WAY 5Y.KO Manual Mass Comparator
UYA 5Y Ultra-Microbalance
MA 5Y.IC Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator
AKM-2 Automatic Mass Comparator
CY10 Precision balance
Analytical Balance for Filters XA 5Y.F

Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
HRP 5Y.KB Mass Comparator
XA 5Y.M Microbalance
RMC F Robotic Weighing System
RB F Robotic Weighing System
UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance for stents weighing
AK-6 5Y.F Automatic Weighing System