



More information on the website  
[radwag.com/us/info,w1,VML](https://radwag.com/us/info,w1,VML)

## 500 $\mu$ l ÷ 5000 $\mu$ l - RW8-107-20-9 Automatic Pipettes



The drawings, photos and graphics used are for illustrative purposes only.

### Datasheet

Volume	500 $\mu$ l ÷ 5000 $\mu$ l
Systematic error	2,0% / 0,5%
Random error	0,6% / 0,15%
Tips	5000 $\mu$ l



### Compatible with

XA 5Y.M.A.P Microbalance  
Microbalance for Pipette Calibration MYA 5Y.P

Automatic Balances for. Pipettes Calibration