



More information on the website
radwag.com/en/info,w1,W13

RES Radwag Extended SI Module (X2)

WX-016-0265



The drawings, photos and graphics used are for illustrative purposes only.

Overview

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MA X2.IC.A Moisture Analyzers
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balances

PS X7 Balances for PGC
PS X2 Balances for PGC
MA X2.A Moisture Analyzers
MA X7.A Moisture Analyzer
PS X7 Precision Balance