




More information on the website
radwag.com/es/info,w1,W13

Protocolo de comunicación Radwag (X2)

WX-016-0265



The drawings, photos and graphics used are for illustrative purposes only.

Descripción

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

Trabaja con (Additional Fee)

Analizador de Humedad
MA X7.IC.A
Balanza analítica AS X7
Balanzas PS X7 para CCE

Analizador de Humedad
MA X7.A
Balanza de precisión PS X7