



More information on the website
radwag.com/en/info,w1,2FA

MAS 220.Y Weighing Module



MAS 220.Y Weighing Module

The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

	MAS 220.Y Weighing Module
Metrological parameters	
Maximum capacity [Max]	220 g
Readability [d]	0,1 mg
Tare range	-220 g
Physical parameters	
Display	5,7" Resistiv-Farb-Touchscreen
Weighing pan dimensions	ø42 mm
Device dimensions	175×236×105 mm
Packaging dimensions	610×340×295 mm
Net weight	3,7 kg
Communication interface	
Communication interface	2×RS232, 2×USB-A, Ethernet, 4 IN / 4 OUT (digital), Wi-Fi
Electrical parameters	
Power supply	100 – 240 V AC 50/60 Hz
Environmental conditions	
Operating temperature	+10 ÷ +40 °C
Relative humidity	40% ÷ 80%

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.



Accessories

Antivibration Tables
Professional Weighing Tables
RS 232, RS 485 cables

RS 232 – USB Converter
RS 232 cables (scale - printer)

Software

RAD-KEY
R-LAB

R Panel
RADWAG Development Studio

Device dimensions

MAS 220.Y Weighing Module

