



SAP/C Industrial Antivibration Table, SAP/H Stainless Steel Industrial Antivibration Table

More information on the website  
radwag.com/en/info,w1,9TS



SAP/C Industrial Antivibration Table



SAP/H Stainless Steel Industrial Antivibration Table

The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

	SAP/C Industrial Antivibration Table	SAP/H Stainless Steel Industrial Antivibration Table
<b>Construction</b>		
<b>Construction</b>	mild steel St3S	Stainless steel AISI304
<b>Load-bearing capacity</b>	40 kg	40 kg
<b>Physical parameters</b>		
<b>Packaging dimensions</b>	820×681×816 mm	820×681×816 mm
<b>Net weight</b>	55 kg	55 kg
<b>Gross weight</b>	72 kg	72,5 kg
<b>Length</b>	770 mm	770 mm
<b>Width</b>	650 mm	650 mm
<b>Height</b>	815±10 mm	815±10 mm


## Compatible with

H315 HR Waterproof Platform Scale  
H315 H Waterproof Platform Scale  
HY10.HRP.H Stainless steel High Resolution Scale  
Stainless Steel HRP EX.H Platforms  
HRP Platforms High Capacity  
HY10.HRP High Resolution Scale

Stainless steel HRP.H platforms  
WPY Multifunctional Scale  
WPY Scales for PGC  
HRP EX Platforms High Capacity  
PUE 7.1.HRP High Resolution Scale  
WLY Scales for PGC

## Device dimensions


### SAP/C Industrial Antivibration Table



**SAP/C**

---

### SAP/H Stainless Steel Industrial Antivibration Table



**SAP/H**

---