

# Protective Covers for Balances

Protection against scratches and dirt.



**A1 protective cover**  
for PS and WLC precision balances



**A2 protective cover**  
for PS and WLC precision balances



**A2 protective cover**  
for PS R.M precision balances



**X2 protective cover**  
for X2 balances

## Features

### Protection

Protective covers secure balances against scratches and dirt. They are particularly useful in harsh ambient conditions under which the balance is exposed to different forms of contamination. The covers are made of transparent plastic and are reusable.

### Various Models

Various models of the covers enable protecting either the display or the whole housing.

## Compatible with

### A1 protective cover for PS and WLC balances

- PS R, WLC A1 precision balances (128×128 mm weighing pan)

### A2 protective cover for PS and WLC balances

- PS R, WLC A2 precision balances (195×195 mm weighing pan)

### A2 protective cover for PS R2.M balance

- PS R.M precision balances

### X2 protective cover:

- AS X2 analytical balances
- PS X2, WLC X2 precision balances
- MA X2 moisture analyzers