

# Voluminous sample container

The set is an ideal solution for drying large-volume samples, e.g. textiles, or products that placed on the weighing pan tend to increase their volume. The container allows to enclose the sample within a limited space wherein the possibility to determine the moisture content value remains unaffected.

*The accessory gives the user full control over the sample and an error-free result.*



# Voluminous sample container



In the case of large-volume samples, the set prevents the sample from coming into contact with solid components of the moisture analyzer (such as the heating element or its cover) and protects it from damage.

Closing a large-volume, low-density sample in the container also reduces the risk of errors resulting from the effect of the sample being raised by warm air. The set consists of 2 parts made of steel.



**MA 3Y**

Max: 60 g - 200 g  
d: 0.1 g - 1 mg



**MA X2**

Max: 50 g - 210 g  
d: 0.1 mg - 1 mg



**MA R**

Max: 50 g - 210 g  
d: 0.1 mg - 1 mg