



MICRO-KIT


MICRO-KIT

Juego de soportes para microtubos

Los soportes especiales de las balanzas RADWAG permiten pesar una amplia gama de microtubos. Están diseñados para aumentar la ergonomía y minimizar el riesgo de contaminación durante la aplicación de la muestra.




Los soportes tienen la posibilidad de inclinar el microtubo, lo que permite colocarlo en cualquier ángulo, dependiendo de la aplicación. El juego incluye corta aires especiales, que estabilizan las condiciones de medición en la cámara de pesaje.



Ejemplos de aplicación

soportes para microtubos



Los diseñadores no olvidaron del organizador de los soportes, lo que permite un fácil almacenamiento de soportes no utilizados.





SENSOR CAMERA SENSOR

Weighing
Fast Dosing

John Smith

Wed, 10.02.2022 09:31:15

4.895774 g

BWP

Net - I - 0 Max

Net weight 4.895774 g
Gross weight 6.534318 g
Tare 11.430092 g
Product name Liquid
Minimum sample weight 0.003 g
MSW status Good

Temperature 19°C

Air Humidity 45%
Air Pressure 980.1 hPa
Air Density 1.158 kg/m3

Vibrations 13%

Navigation icons: Settings, Home, 123, Bar chart, Tare, Tare, Tare, Settings, Tare, Home

Home >0< >T< Home



www.radwag.com