



You Are Safe Now

And so is your terminal.

The screen of the 5Y and CY10 terminal is made of glass – a material prone to scratches and cracks. It can easily be dropped or damaged by transport, so it is worth ensuring that it is properly protected.

Safety in Every Inch

Hybrid glass protects the touch panel of the terminal from scratches and cracks. It combines the advantages of glass and film: resistance and flexibility. It is perfectly matched to the balance's or scale's display, so it provides the highest level of security to every inch of the screen.

Instant Responsiveness

The hybrid protective glass allows instant response to touch. It does not affect the user experience of the terminal – it does not limit the sensitivity of the screen, does not reduce viewing angles and does not distort colors.

Hybrid Protective Glass at a Glance:

- very thin – ~0.3 mm
- perfectly smooth
- durable – strengthens the screen by up to 250%
- scratch-resistant – glass hardness of 7H

The Set Includes:

- hybrid protective glass
- leave that removes pollen from under the already applied glass
- positioning sticker to help fit the glass to the screen
- cardboard to eliminate air bubbles
- cleaning kit (moistened cloth and non-dusting cloth)

You can buy and install the hybrid protective glass at any time. Installation is very easy and takes a while.



Hybrid of Resistance & Flexibility
Hybrid glass combines the advantages of glass and film: it will neither scratch nor chip.



The Peak of Security
Strengthening the screen by 250% means that the terminal will emerge from a fall from a height to a hard surface – unscathed.



So Thin, yet so Strong
This almost invisible and undetectable coating translates not only into solid protection, but also into a pleasure to use.



Off-Scale Protection Level
Glass hardness of 7H guarantees you a screen without a single scratch.

