

# Ethernet Cables

Compatibility of weighing instruments and indicators, and integration with peripheral equipment.



P0198 cable



PT0212 cable



PT0303 cable



PT0302 cable

## Compatible With

### P0198

- Weighing instruments with PUE C41H indicator
- Weighing instruments with PUE HY indicator
- PUE C41H indicator
- PUE HY indicator
- HRP platforms
- MWSH weighing module
- MWLH weighing module

### PT0212

- Weighing instruments with PUE 5.15 and 5.19 indicator
- PUE 5.15 and 5.19 indicator

### PT0302

- HRP platforms
- MWSH weighing module
- MWLH weighing module

### PT0303

- HRP platforms
- MWSH weighing module
- MWLH weighing module

## Technical Specifications

	P0198	PT0212	PT0302	PT0303
<b>Connector</b>	M12 4P	RJ45	weighing instrument (HRP, MWSH, MWLH) – indicator (HY10, PUE5)	weighing instrument (HRP, MWSH, MWLH) – indicator (HY10, PUE5)
<b>Length</b>	3 m	3 m	5 m	5 m

## Dimensions



### P0198