

Protective Covers for Balances

Protection against scratches and dirt.



A1 protective cover for PS and WLC precision balances with 128x128 mm weighing pan



A2 protective cover for PS and WLC balances with 195 x 195 mm weighing pan



A2 protective cover for PS R2.M precision balance



X2 protective cover

Features

Protection

Protective covers secure balances against scratches and dirt. They are particularly useful in harsh ambient conditions under which the balance is exposed to different forms of contamination. The covers are made of transparent plastic and are reusable.

Various Models

Various models of the covers enable protecting either the display or the whole housing.

Compatible With

A1 protective cover for PS and WLC balances with 128x128 mm weighing pan:

- PS R, WLC precision balances

A2 protective cover for PS and WLC balances with 195 x 195 mm weighing pan:

- Wag precyzyjnych PS R; WLC, C/2

A2 protective cover for PS R2.M precision balance

- PS R precision balance

X2 protective cover:

- AS X2 analytical balances
- PS X2, WLC precision balances