

MICRO-KIT

Set of holders for microscale glassware

Special holders available in RADWAG balances enable weighing of a wide range of microscale glassware. They have been designed to increase ergonomics and minimize the risk of contamination during sample application.



The holders have the ability to tilt the microscale glassware, allowing it to be placed at any angle, depending on the application. The set includes dedicated draft shields that stabilize the measurement conditions in the weighing chamber.





SET OF HOLDERS FOR MICROSCALE GLASSWARE

Examples of application of the microscale glassware holders













