



# IM02.1 (MediaBox) Communication Module

WX-016-0303

More information on the website  
[radwag.com/en/info,w1,07X](http://radwag.com/en/info,w1,07X)



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Communication interface	
Communication interface	RS232, 4IN/4OUT, Virtual Com
Physical parameters	
Device dimensions	100×100×30.5 mm




## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 5Y.PM Precision Balance  
 MYA 5Y.F Microbalance for Filters  
 XA 5Y.M.A Microbalance  
 CY10 Scales for PGC  
 PUE CY10 Weighing Terminal  
 CY10 Multifunctional Scale  
 XA 5Y.A Analytical Balance  
 UYA 5Y Ultra-Microbalance

CY10 Precision Balance  
 XA 5Y.F Analytical Balance for Filters  
 MYA 5Y Microbalance  
 XA 5Y Analytical Balance  
 PMV 5Y Moisture Analyzer  
 XA 5Y.M Microbalance  
 RS 232, RS 485 cables  
 UYA 5Y.F Ultra-Microbalance for Filters  
 XA 5Y.M.A.S Microbalance for stents weighing

## Device dimensions



100×100×30,5 mm