



More information on the website
radwag.com/en/info,w1,080

AK-5/100.STENT Automatic Balance for Stents



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------------------|--------------------------------|
| Maximum capacity [Max] | 5.1 g |
| Readability [d] | 0.1 µg |
| Tare range | -5.1 g |
| Standard repeatability [5% Max] | 0.8 µg |
| Standard repeatability [Max] | 1 µg |
| Standard minimum weight (USP) | 0.8 mg |
| Standard minimum weight (U=1%, k=2) | 0.08 mg |
| Permissible repeatability [5% Max] | 1.5 µg |
| Permissible repeatability [Max] | 2 µg |
| Linearity | ±5 µg |
| Adjustment | external |
| Physical parameters | |
| Display | 10" graphic colour touchscreen |
| Device dimensions | 395×195×420 mm |
| Packaging dimensions | 860×800×550 mm |

| Physical parameters | |
|-------------------------|--|
| Net weight | 20 kg |
| Gross weight | 30 kg |
| Communication interface | |
| Communication interface | 2×USB-A, USB-C, HDMI, Ethernet, Wi-Fi, Hotspot |
| Electrical parameters | |
| Power supply | 100 – 240 V AC 50/60 Hz |

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.



Accessories

Receipt Printer
Fingerprint Reader

Additional modules