



More information on the website
radwag.com/en/info,w1,0AD

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 10 kg), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 21 110 g |
| Nominal value | 1 g – 10 kg |
| Quantity | 17 pcs |

Compatible with

Tweezers

Brushes