




More information on the website  
[radwag.com/en/info,w1,1KQ](http://radwag.com/en/info,w1,1KQ)

# RAD Key



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Components and software

Operating system	Windows 7, 8, 10
------------------	------------------

## Compatible with

H315 HR Waterproof Platform Scale  
 MA X7.IC.A Moisture Analyzer  
 MUYA Modules  
 HX7 H Stainless steel Multifunctional Scale  
 MAS1 Modules  
 HX7.4N H Multifunctional Stainless Steel Ramp Scale  
 HX7.4C Multifunctional Platform Scale  
 HX7 HR Multifunctional Scale  
 XA 5Y.M.A.P Microbalance for pipettes calibration  
 WPY Multifunctional Scale  
 Microbalance for Pipette Calibration MYA 5Y.P  
 C32 Multifunctional Scale  
 Stainless Steel HRP EX.H Platforms  
 HX5.EX-1.4P2 H Stainless Steel Beam Scales  
 HY10.HRP.H Stainless steel High Resolution Scale  
 MA 5Y Moisture Analyzer  
 C315 Platform scale  
 Microbalance for Filters MYA 5Y.F

H315 4N.H.LD Stainless Steel Ramp Scale  
 PS X7 Precision Balance  
 WPY Scales for PGC  
 XA 5Y.A Analytical Balance  
 HX7.4 H Multifunctional Stainless Steel Platform Scale  
 C315.4P pallet scales  
 HX7.4P2 Multifunctional Beam Scales  
 WAY 5Y.KO Manual Mass Comparator  
 UYA 5Y Ultra-Microbalance  
 MA 5Y.IC Moisture Analyzer  
 HX5.EX-1 C 1-Load-Cell Platform Scales  
 MA R Moisture Analyzer  
 HX5.EX-1 H Waterproof Platform Scales  
 AS 5Y Analytical Balance  
 AK-4 Automatic Mass Comparator  
 AKM-2 Automatic Mass Comparator  
 PS R2.H Precision Balance  
 PUE HX5.EX Weighing Terminal

HX5.EX-1 HR Stainless Steel 1-Load-Cell Platform Scales  
HX7.4P2 H Multifunctional Stainless Steel Beam Scales  
TMX H Stainless steel Multifunctional Scale  
WLY Scales for PGC  
MA X2.IC.A Moisture Analyzer  
AS X7 Analytical Balance  
PS X2 Precision Balance  
HRP Platforms High Capacity  
AS X2 PLUS Analytical Balance  
TMX Multifunctional Scale  
C315.D Baby Scale  
MWMH Modules  
XA 5Y.M.A Microbalance  
WLC Precision Balance  
PS R2 Precision Balance  
C315.4P2 Beam Scales  
H315 4H Stainless Steel Platform Scale  
H315 4P2.H Stainless steel Beam Scales  
WLC X7 Precision Balance  
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale  
HX5.EX-1.4P2 C Beam Scales  
PS X2 Balances for PGC  
Stainless steel HRP.H platforms  
H315 H Waterproof Platform Scale  
PUE CY10 Weighing Terminal  
CY10 Multifunctional Scale  
HX7.4P H Multifunctional Stainless Steel Pallet Scales  
MA X2.A Moisture Analyzer  
HX5.EX-1.4N Stainless Steel Ramp Scales  
HY10 Multifunctional Scale  
HX5.EX-1.4 C 4-Load-Cell Platform Scales  
MAS Modules  
HY10 H Stainless Steel Multifunctional Scale  
H315.4 H/Z Stainless Steel Platform Scale, pit version  
MPS Modules  
H315 4N.H Stainless Steel Ramp Scale  
HRP EX Platforms  
C315.4 Platform Scale  
HX5.EX-1.4P H Stainless Steel Pallet Scales  
HY10.HRP High Resolution Scale  
C32.PM Precision Balance  
PUE 7.1.HRP High Resolution Scale  
Automatic Balances for. Pipettes Calibration  
HX7 Multifunctional Scale  
MWLH Modules  
HX5.EX-1.4 H/Z Stainless Steel Platform Scales, Pit Version  
MA X7.A Moisture Analyzer  
HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version  
C315.O Personal Scale

HY10.H Scales for PGC  
HX5.EX-1.4 H Stainless Steel 4-Load-Cell Platform Scales  
WTC Precision Balance  
CY10 Precision Balance  
H315 4C Platform Scale  
C315.OW Personal Scale  
XA 5Y.F Analytical Balance for Filters  
WLY Precision Balance  
UMA Automatic Mass Comparator  
H315 2K Overhead Track Scales  
APP 5Y.KO Manual Mass Comparator  
H315 4P.H Stainless Steel Pallet Scales  
Robotic Mass Comparators  
MW-01 Mass Converters  
WLC C/2 Precision Balance  
MYA 5Y Microbalance  
H315 4I.S Livestock Scales  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
MA X2.A Moisture Analyzers  
5Y.PM.KB Manual Mass Comparator  
AS R2 PLUS Analytical Balance  
PUE C32 Weighing Terminal  
PUE 7.1 Weighing Terminal  
PUE HX7 Weighing Terminal  
PUE 5 Weighing Terminal  
HRP 5Y.KO Mass Comparator  
PS 5Y Precision Balance  
Postal Scales for Packages  
PMV 5Y Moisture Analyzer  
HX5.EX-1.4P C Pallet Scales  
XA 5Y.KO Manual Mass Comparator  
HRP 5Y.KB Mass Comparator  
H315 4I.H Stainless Steel Livestock Scales  
XA 5Y.M Microbalance  
PUE C315 Weighing Terminal  
PUE HY10 Weighing Terminal  
MW-04 Mass Converters  
UYA 5Y.KO Manual Mass Comparator  
AS R1 PLUS Analytical Balance  
Ultra-Microbalance for Filters UYA 5Y.F  
C315.4B Bed Scale  
PUE H315 Weighing Terminal  
HX7.4P Multifunctional pallet scales  
MWSH Modules  
Nano-Comparator NANO.AK-4.500.5Y  
XA 5Y.M.A.S Microbalance for stents weighing  
H315 Stainless Steel Overhead Track Scales  
C315.8B Ramps Bed Scale