



More information on the website
radwag.com/tr/info,w1,10F

UL-04-1 (AS) Tüp Tutucu



The drawings, photos and graphics used are for illustrative purposes only.

Teknik Veriler

Fiziksel parametreler

| | |
|-------------------|---------------|
| Ambalaj boyutları | 110×102×56 mm |
| Net ağırlık | 103 g |
| Brüt ağırlık | 145 g |

İle uyumlu

AS X7 Analitik Teraziler
AS X2 PLUS Analitik Teraziler
AS 5Y Analitik Teraziler

AS R2 PLUS Analitik Teraziler
AS R1 PLUS Analitik Teraziler