



More information on the website  
[radwag.com/en/info,w1,2HV](https://radwag.com/en/info,w1,2HV)

# Knob weights, Set (1 g - 200 g), Wooden Box

OK-101-0032



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class E2

### Construction

Material Stainless steel

Weight 610 g

Nominal value 1 – 200 g

Quantity 11 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes