



More information on the website
radwag.com/en/info,w1,3KF

ZR-04/1 Powerbank



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Electrical parameters | |
|-----------------------|--|
| Power supply | Type-C: 5V/3A, 9V/3A, 12V/3A, 15V/3A, 20V/3A 60W (max) (ZR-04) |
| Physical parameters | |
| Device dimensions | 153×69×42 mm |
| Packaging dimensions | 230×190×80 mm |
| Net weight | 645 g |
| Gross weight | 900 g |

Model: Baseus Ambilight 65W Power bank 30000 mAh



Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 C32 Multifunctional Scale

C315.4P pallet scales
 UYA 5Y Ultra-Microbalance
 AS 5Y Analytical Balance

C315 Platform scale
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
WLY Scales for PGC
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balance
PS X7 Balances for PGC
C315.D Baby Scale
XA 5Y.M.A Microbalance
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
WLC X7 Precision Balance
CY10 Scales for PGC
PS X2 Balances for PGC
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
C315.4 Platform Scale
C32.PM Precision Balance
PS X7 Precision Balance
WPY Scales for PGC
XA 5Y.A Analytical Balance

PS 5Y Balances for PGC
PS R2.H Precision Balance
WTC Precision Balance
CY10 Precision Balance
C315.OW Personal Scale
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
WLC C/2 Precision Balance
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PS 5Y Precision Balance
XA 5Y.M Microbalance
PUE C315 Weighing Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed Scale
C315.K Chair Scale
XA 5Y.M.A.S Microbalance for stents weighing
C315.8B Ramps Bed Scale