



More information on the website  
[radwag.com/en/info,w1,3Y5](https://radwag.com/en/info,w1,3Y5)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 1 kg), Plastic Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	2 111.11 g
Nominal value	1 mg - 1 kg
Quantity	25 pcs

## Compatible with

Tweezers

Brushes