



More information on the website  
[radwag.com/en/info,w1,5HY](http://radwag.com/en/info,w1,5HY)

# Knob weights, Set (1 g - 5 kg), Wooden Box

OK-101-0036



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	11 110 g
Nominal value	1 g – 5 kg
Quantity	16 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes