



More information on the website  
[radwag.com/en/info,w1,8BD](https://radwag.com/en/info,w1,8BD)

# F2 Mass Standard - Knob Weights - 2000 kg - Service Case

OP-304-0012



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F2              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 2000 kg         |

## Compatible with (Additional Fee)

Brushes