



More information on the website
radwag.com/en/info,w1,8PJ

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 10 kg), Plastic Box

OK-302-0041



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	21 111.11 g
Nominal value	1 mg – 10 kg
Quantity	29 pcs

Compatible with (Additional Fee)

Tweezers

Brushes