




More information on the website
radwag.com/it/info,w1,8RT

Protocollo di comunicazione Radwag (5Y, CY10)



The drawings, photos and graphics used are for illustrative purposes only.

Caratteristiche generali

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

Compatibile con

Microbilance XA 5Y.M.A.P
Microbilancia per la calibrazione delle pipette MYA 5Y.P
Sistema di pesatura automatico UMA 5Y.F
Termobilance MA 5Y
Bilance di precisione 5Y.PM
Microbilancia per filtri MYA 5Y.F
Microbilance XA 5Y.M.A

Comparatore di massa automatico UMA
APP 5Y.KO Comparatore di massa manuale
Comparatori di densità
Microbilance MYA 5Y
Bilance analitiche XA 5Y
5Y.PM.KB Comparatore di massa manuale
HRP 5Y.KO Comparatore di massa manuale

Terminale di pesata PUE CY10
Bilance multifunzione CY10
Susceptometro per la misurazione del magnetismo
Bilance automatiche per la taratura delle pipette
Bilance analitiche XA 5Y.A
WAY 5Y.KO Comparatore di massa manuale
Comparatore di massa robotici
UYA 5Y Ultra-Microbilance
Termobilance MA 5Y.IC
Bilance analitiche AS 5Y
Comparatore di massa automatico AK-4
Comparatore di massa automatico AKM-2
Bilance di precisione CY10
Bilancia analitica per filtri XA 5Y.F

Bilance di precisione PS 5Y
Comparatore di massa automatico
Termobilance PMV 5Y
XA 5Y.KO Comparatore di massa manuale
HRP 5Y.KB Comparatore di massa manuale
Microbilance XA 5Y.M
Sistema di pesatura robotizzato RMC F
Sistema di pesatura robotizzato RB F
UYA 5Y.KO Comparatore di massa manuale
Ultra-microbilance per filtri UYA 5Y.F
Nano-Comparatore NANO.AK-4.500.5Y
Microbilance XA 5Y.M.A.S
Sistema di pesatura automatico AK-6 5Y.F