



More information on the website  
[radwag.com/en/info,w1,8V1](http://radwag.com/en/info,w1,8V1)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 500 g), Plastic Box

OK-202-0057



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 1 110 g         |
| Nominal value           | 1 – 500 g       |
| Quantity                | 12 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes