



More information on the website
radwag.com/de/info,w1,9KC

Knopfgewichte, Set (1 mg - 200 g), Holz-Etui

OK-401-0081



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter

| | |
|-------------|--------------|
| OIML-Klasse | M1 |
| Gewicht | 611,11 g |
| Nominalwert | 1 mg – 200 g |
| Anzahl | 23 St. |

Kompatibel mit (Additional Fee)

Pinsel