



More information on the website  
[radwag.com/en/info,w1,BJ7](https://radwag.com/en/info,w1,BJ7)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 5 g), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 6.11 g          |
| Nominal value           | 1 mg - 5 g      |
| Quantity                | 16 pcs          |

## Compatible with

Tweezers

Brushes