



More information on the website
radwag.com/us/info,w1,BS0

Hybrid Protective Glass 5Y, CY10

WX-001-0435



The drawings, photos and graphics used are for illustrative purposes only.

Overview

Protective cover is made of transparent plastic, secures a display against scratches, dust etc. It is intended for balances X2 and R2 type.

Set - 10 pcs.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 UMA 5Y.F Automatic Weighing System
 5Y.PM Precision balance
 MYA 5Y.F Microbalance for Filters
 XA 5Y.M.A Microbalance
 CY10 Scales for PGC
 PUE CY10 Weighing Terminal
 CY10 Multifunctional Scale
 Automatic Balances for. Pipettes Calibration
 XA 5Y.A Analytical Balance
 WAY 5Y.KO Manual Mass Comparator
 UYA 5Y Ultra-Microbalance
 AS 5Y Analytical Balance
 AK-4 Automatic Mass Comparator
 PS 5Y Balances for PGC

Robotic Mass Comparators
 Density measurement
 MYA 5Y Microbalance
 XA 5Y Analytical Balance
 Manual susceptometers - Magnetization
 5Y.PM.KB Manual Mass Comparator
 HRP 5Y.KO Mass Comparator
 PS 5Y Precision Balance
 Automatic Vacuum Mass Comparators
 PMV 5Y Moisture Analyzer
 Automatic susceptometers - Magnetization
 XA 5Y.KO Manual Mass Comparator
 HRP 5Y.KB Mass Comparator
 XA 5Y.M Microbalance
 RMC F Robotic Weighing System
 RB F Robotic Weighing System

AKM-2 Automatic Mass Comparator
CY10 Precision balance
Analytical Balance for Filters XA 5Y.F
UMA Automatic Mass Comparator
APP 5Y.KO Manual Mass Comparator

UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance for stents weighing
AK-6 5Y.F Automatic Weighing System