



More information on the website
radwag.com/us/info,w1,BS0

Hybrid Protective Glass 5Y, CY10



The drawings, photos and graphics used are for illustrative purposes only.

Overview

Protective cover is made of transparent plastic, secures a display against scratches, dust etc. It is intended for balances X2 and R2 type.

Set - 10 pcs.

Compatible with

XA 5Y.M.A.P Microbalance
 Microbalance for Pipette Calibration MYA 5Y.P
 UMA 5Y.F Automatic Weighing System
 5Y.PM Precision balance
 Microbalance for Filters MYA 5Y.F
 XA 5Y.M.A Microbalance
 PUE CY10 Weighing Terminal
 CY10 Multifunctional Scale
 Susceptometer - Magnetization
 Automatic Balances for Pipettes Calibration
 XA 5Y.A Analytical Balance
 WAY 5Y.KO Manual Mass Comparator
 UYA 5Y Ultra-Microbalance
 AS 5Y Analytical Balance
 AK-4 Automatic Mass Comparator
 AKM-2 Automatic Mass Comparator

Robotic Mass Comparators
 Density measurement
 MYA 5Y Microbalance
 XA 5Y Analytical Balance
 5Y.PM.KB Manual Mass Comparator
 HRP 5Y.KO Mass Comparator
 PS 5Y Precision Balance
 Automatic Vacuum Mass Comparators
 PMV 5Y Moisture Analyzer
 XA 5Y.KO Manual Mass Comparator
 HRP 5Y.KB Mass Comparator
 XA 5Y.M Microbalance
 RMC F Robotic Weighing System
 RB F Robotic Weighing System
 UYA 5Y.KO Manual Mass Comparator
 Ultra-Microbalance for Filters UYA 5Y.F

CY10 Precision balance
Analytical Balance for Filters XA 5Y.F
UMA Automatic Mass Comparator
APP 5Y.KO Manual Mass Comparator

Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance
AK-6 5Y.F Automatic Weighing System