



More information on the website
radwag.com/en/info,w1,BSB

LS2208 Barcode Scanner



The drawings, photos and graphics used are for illustrative purposes only.

Overview

LS2208 barcode scanner is a manual device designed for operation in adverse industrial conditions. It integrates with RADWAG terminals, PUE H315, PUE 5, PUE 7.1, HY10, PUE C32, and balances, 5Y, R2, X2, H7.

LS2208 is a high quality handheld barcode scanner. It has been designed for intense and continuous operation in retail sales, libraries, health centers and pharmacy, etc. Its light and ergonomic casing provides comfort of use. The function of dynamic scanning and wide scanning range, enable LS2208 to collect barcodes in a quick and reliable way, increasing work pace and efficiency. Intellistand™ system provides scanning of codes placed in front of the scanner, and easy switching between handheld and handfree mode (operation on a stand). A patented design of a single printed-circuit board eliminates failure endangered parts, i.e. connections between multiple boards. Scanning part of LS2208 with liquid polymere protects its module from shock based damages. LS2208 scanner features like integrated operation of various systems, universal cable or new symbols detection guarantee indispensable protection of an investment now and in the future.

Compatible with

MA X7.IC.A Moisture Analyzer
MUYA Modules

XA 5Y.A Analytical Balance
HX7.4 H Multifunctional Stainless Steel Platform Scale

HX7 H Stainless steel Multifunctional Scale
 HX7.4N H Multifunctional Stainless Steel Ramp Scale
 HX7.4C Multifunctional Platform Scale
 HX7 HR Multifunctional Scale
 XA 5Y.M.A.P Microbalance for pipettes calibration
 WPY Multifunctional Scale
 MYA 5Y.P Microbalance for pipettes calibration
 C32 Multifunctional Scale
 HY10.HRP.H Stainless steel High Resolution Scale
 MA 5Y Moisture Analyzer
 5Y.PM Precision Balance
 MYA 5Y.F Microbalance for Filters
 HX7.4P2 H Multifunctional Stainless Steel Beam Scales
 TMX H Stainless steel Multifunctional Scale
 WLY Scales for PGC
 MA X2.IC.A Moisture Analyzer
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balance
 PS X7 Balances for PGC
 TMX Multifunctional Scale
 XA 5Y.M.A Microbalance
 PS R2 Precision Balance
 WLC X7 Precision Balance
 HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale
 CY10 Scales for PGC
 PS X2 Balances for PGC
 DWT/RC/HYF Checkweighers
 PUE CY10 Weighing Terminal
 CY10 Multifunctional Scale
 HX7.4P H Multifunctional Stainless Steel Pallet Scales
 MA X2.A Moisture Analyzer
 HY10 Multifunctional Scale
 HY10 H Stainless Steel Multifunctional Scale
 HY10.HRP High Resolution Scale
 C32.PM Precision Balance
 PUE 7.1.HRP High Resolution Scale
 Automatic Balances for. Pipettes Calibration
 HX7 Multifunctional Scale
 MA X7.A Moisture Analyzer
 HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version
 PS X7 Precision Balance
 WPY Scales for PGC
 HX7.4P2 Multifunctional Beam Scales
 WAY 5Y.KO Manual Mass Comparator
 UYA 5Y Ultra-Microbalance
 MA 5Y.IC Moisture Analyzer
 MA R Moisture Analyzer
 AS 5Y Analytical Balance
 AK-4 Automatic Mass Comparator
 PS 5Y Balances for PGC
 AKM-2 Automatic Mass Comparator
 PS R2.H Precision Balance
 HY10.H Scales for PGC
 CY10 Precision Balance
 XA 5Y.F Analytical Balance for Filters
 WLY Precision Balance
 UMA Automatic Mass Comparator
 APP 5Y.KO Manual Mass Comparator
 Robotic Mass Comparators
 Density Measurement
 MYA 5Y Microbalance
 PS R1 Precision Balance
 XA 5Y Analytical Balance
 MA X2.A Moisture Analyzers
 5Y.PM.KB Manual Mass Comparator
 AS R2 PLUS Analytical Balance
 PUE C32 Weighing Terminal
 PUE 7.1 Weighing Terminal
 PUE HX7 Weighing Terminal
 PUE 5 Weighing Terminal
 HRP 5Y.KO Mass Comparator
 PS 5Y Precision Balance
 Automatic Vacuum Mass Comparators
 PMV 5Y Moisture Analyzer
 XA 5Y.KO Manual Mass Comparator
 HRP 5Y.KB Mass Comparator
 XA 5Y.M Microbalance
 UYA 5Y.KO Manual Mass Comparator
 AS R1 PLUS Analytical Balance
 UYA 5Y.F Ultra-Microbalance for Filters
 H315 H Waterproof Platform Scale
 HX7.4P Multifunctional pallet scales
 Nano-Comparator NANO.AK-4.500.5Y
 XA 5Y.M.A.S Microbalance for stents weighing