



More information on the website
radwag.com/en/info,w1,CBG

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 100 g), Wooden Box

OK-201-0054



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 211.11 g |
| Nominal value | 1 mg – 100 g |
| Quantity | 21 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes