



Transport Case for PM 4Y.KB Mass Comparator



Overview

The case has been designed for safe transportation of the PM 4Y.KB mass comparator. It is equipped with two handles for convenient carrying so that the device can be transported by two people at the same time (recommended). The interior of the case has been adapted to fit all device components, with this it can be carried without the risk of damage.

The case is made of 6 mm thick plywood. It features reinforced corners, aluminium trim, and hard dark foam filling.

The case comes as an accessory.

Compatible with

PM.4Y.KB Manual Mass Comparator