



More information on the website
radwag.com/en/info,w1,CXS

SA/APP/C - Antivibration Table for Mass Comparators



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Construction | |
|-----------------------|-----------------|
| Construction | mild steel St3S |
| Load-bearing capacity | 150 kg |
| Physical parameters | |
| Packaging dimensions | 710×470×225 mm |
| Net weight | 128 kg |
| Gross weight | 142 kg |
| Length | 710 mm |
| Width | 470 mm |
| Height | 225 mm |

Compatible with

APP 5Y.K0 Manual Mass Comparator

5Y.PM.KB Manual Mass Comparator

Device dimensions

