



More information on the website
radwag.com/en/info,w1,DXV

Transport Case for MA X2, MAR Moisture Analyzer

WX-003-0023



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters

Packaging dimensions	545×330×340 mm
----------------------	----------------

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MA 5Y Moisture Analyzer
MA X2.IC.A Moisture Analyzers
MA X2.A Moisture Analyzers

MA X7.A Moisture Analyzer
MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer

Device dimensions

